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Surgical instrumentation must be used in premises compliant with standard aseptic conditions and with current practices in health care establishments The ancillaries and instruments are delivered clean but NON-STERILE. The transfer of an ancillary between two health care establishments is forbidden without prior return to ADSM.

Certain ancillaries are constituted by temporary invasive devices such as drills, reamers, mills, taps and probes. These instruments must be used with precaution. The equipment must only be used for the function for which they are dedicated in operating techniques. Any deterioration of the equipment can lead to risks of malfunction

GENERALITIES

Before being used, the instruments, whether new or not, must be treated carefully according to the instructions. The instruments boxes do not constitute a sterile barrier; a sterile packaging must be used to keep them in a sterile state Ancillary equipment is delivered non-sterile in transport containers. A transfer sheet (which, according to the final destination before being sent back to ADSM) is supplied with each ancillary for checking at reception and then before sterilisation

Instruments or ancillaries generally have a long service life. However, poor handling or unsuitable production methods reduce the life of the instrument, which may be the cause of long use, poor handling or inappropriate maintenance, must be returned to ADSM. Any malfunction concerning the instruments must be indicated to ADSM ADSM instruments and ancillaries must undergo all the steps of decontamination, cleaning, inspection and final sterilisation before being sent back to ADSM. The documentation on decontamination must be supplied with instruments returned to ADSM. Any non-sterile medical device received must be inspected upon receipt to verify the integrity and the proper functioning of the medical device. Washer-disinfectors used by the health care facility must comply with ISO 15883.

CLEANING

Caution

Vibrating precautions must be observed by hospital personnel who work with contaminated or potentially contaminated medical devices. Pointed or cutting devices must be handled with care Do not use metal brushes or scouring pads during manual cleaning procedures. These items risk damaging the surface and the finish of the instrument. Use nylon brushes with soft bristles and swabs with dimensions adapted to the devices to treat Do not allow contaminated devices to dry before re-treatment so as to facilitate later steps Saline solution and cleaning/disinfecting agents containing aldehyde, mercury, active chlorine, chloride, bromine, bromide, iodine or iodide are not to be used Do not use mineral oil or silicone lubricants as they coat the micro-organisms preventing direct contact of the steam with the surface and are difficult to eliminate. Cleaning and neutral pH enzymatic agents are recommended and preferred for cleaning reusable ADSM devices. Only agents with a neutral pH must be used for cleaning aluminium devices. Contact with strong alkaline detergents or with solutions containing caustic soda, iodine or chlorine must be avoided as aluminium and stainless steel can be chemically attacked and the device damaged. The manufacturer must always be refered to and comply with the instructions provided by the cleaning product manufacturer

The use of hard water must be avoided. Softened tap water can be used for initial rinsing. Final rinsing must be done with purified water in order to eliminate deposits on the instruments. One of the first steps in the cleaning process can be employed to purify water: ultrafiltration (UF), reverse osmosis (RO), de-ionisation or equivalent Treatment 1: Drills, broaches, rasps and cutting instruments must be carefully inspected after treatment with alkaline detergents in order to ensure that the cutting edges are not deteriorated. Comment 2: In the case of first use of the instruments, the establishment shall perform a series of cleaning operations in order to eliminate any earlier contamination of the instruments. A minimum of 3 consecutive cleaning operations is recommended.

Instructions

When the device will be used Remove excess liquid and organic tissues on the instruments with a lint-free disposable pad. Place the devices in a tray of distilled water and cover with a damp cloth Remove the plastic films and potential protective foam Usually check that all the instruments are present and that the attached check-list are present Comment: the instruments must be cleaned within 30 minutes following their use to minimise the possibility of drying before cleaning. **Confirmation and transport** The instruments used must be sent to the procurement service in closed or covered containers in order to avoid any new risk of contamination When this is applicable, the instruments with several components must be disassembled for efficient cleaning. Make sure not to lose small screws and components When this is applicable, articulated instruments must be kept in their original containers. All cleaning agents must be prepared observing the dilution and the temperature recommended by the manufacturer. Softened tap water can be used for preparing the cleaning agents The use of recommended temperatures is important for optimal performances of the cleaning agents Comment: new cleaning solutions must be prepared.

Procedure for automatic cleaning/disinfection

Soaking (pre-cleaning) Completely immerse the instruments in an alkaline solution having fungicide, bactericides and virucide powers and leave them to soak for at least 15 minutes at room temperature. We recommend the use of the product NEOCHERBERPTO PRECLEAN ZP (C) 5 m/L of the brand DR WEIGERT Use a nylon, flexible bristle brush to softly brush the device until all visible stains are eliminated. Pay particular attention to asperities, cavities, contact surfaces, connectors and other difficult to clean areas. Action the articulated parts in order to get access to all the devices Comment: it is not recommended to prolong the soaking step for an indeterminate period at the risk of deteriorating the instruments. Remove the device from the alkaline solution and rinse abundantly with lukewarm water. The use of hot water for less than 5 minutes. Carefully and effectively rinse cavities, holes and any other difficult access areas

Cleaning Place the instruments in the tray of a suitable washing machine/autoclave and treat using the standard cycle of the instrument washing machine/autoclave per the instructions from the manufacturer of the washing machine/autoclave Washing with a detergent and disinfectant solution suitable for the washing machine/autoclave (observe the contraindications and the temperature of use of the products). We recommend the use of the product NEOCHERBERPTO CLEAN (C)3-5 mL DRG WEIGERT's brand.

Disinfection

Thermal disinfection assuring a value of *t*, greater than 3000, i.e. A minimum of 3 min at 93°C (or any other equivalent time/temperature pair) Rinse with demineralised water Drying not exceeding 110°C. We recommend the use of MEDIKLAR® SPECIAL drying vibrator (C) 8 m of the brand DR WEIGERT's brand Comment: The ancillary equipment is designed to operate at room temperature. It is recommended not to use the equipment directly after steam sterilisation and to wait for the temperature to return to ambient. Before the intervention, it is the responsibility of the person in charge of operating the equipment to check the proper coupling of the ADSM instruments with the equipment available in the operating theatre.

Inspection before packaging

Carefully inspect each device to ensure the equipment offers any visible contamination. If there is contamination, repeat the cleaning operations. Check the operation of the equipment. Check the action of movable parts (e.g. hinges, connectors, sliding parts, springs, etc.) over their complete amplitude Check that there is no deformation of the instruments, in particular instruments which are connected to motors When the instruments are assembled, check that the devices fit together correctly with any components in contact Check that there is no excessive play between the assembled components

PACKAGING

Medical quality packaging of steel sterilisation sachets available commercially, can be used for wrapping individual instruments. The wrapping must be prepared using a protective wrapping technique using an equivalent method For sterilisation, trays and boxes with covers can also be placed in an approved sterilisation container with sealing lid. Follow the sterilisation container manufacturer's instructions for insertion and the replacement of sterilisation filters in sterilisation containers Instrument transport and storage must be done as follows in the following conditions: All devices must be properly arranged in order to ensure contact of the steam with all instrument surfaces. The instruments must not be stacked or placed in close contact (use the arranged supports or silicone mat provided for this) The user must check the contents of the instrument box do not spill out, once the devices are arranged in the box Only devices manufactured and/or distributed by ADSM must be placed in ADSM trays

STERILIZATION

The cleaning/decontamination step is compulsory before sterilising. It is important that sufficient cleaning is done before sterilisation without which, the sterilisation can be compromised After cleaning, the instruments must be rapidly and meticulously reassembled and placed in their respective trays/locations in order to avoid any avoidable damage or contamination The instrument trays and supports must have been decontaminated and cleaned according to the same arrangements as those for the instruments. The health care establishment must be responsible for the internal procedures for reassembly, inspection and packaging of instruments in order to ensure the sterilising penetration of the steam and appropriate drying. The arrangements to make for the protection of sharp or potentially injuring edges of instruments must also be recommended by the health care establishment Sterilisation by steam/most heat is the also recommended by the ADSM set of orthopaedic instruments. Refer to the table below for the parameters of the sterilisation cycles per instrument, making sure to respect the maximum level (SAL) of 10⁻⁶ Sterilisation methods using ethylene oxide or plasma must not be used The recommendations of the manufacturer of the steriliser must always be followed. During the sterilisation of several instrument sets in a sterilisation cycle, ensure that the maximum load indicated by the manufacturer is not exceeded. The set of instruments must be correctly prepared and packaged in the trays and/or boxes to allow the steam to penetrate and enter in direct contact with all surfaces

Validated parameters for steam sterilization:

Type of cycle	Temperature (°C)/(°F)	Exposure Time	Drying Time
US	Pre-vacuum 132°C / 270°F	4 min	20 min
International*	Pre-vacuum 134°C / 273°F	3 min	20 min
World Health Organization (WHO)**	Pre-vacuum 134°C / 273°F	18 min	20 min

NB:

*7/10 steam sterilization cycles with longer times than those listed are also acceptable **/10 This cycle is not recommended in the USA ** Recommended by the for reprocessing instruments where there is concern regarding TSE/CJD contamination. In the case of an operation with a risk of Non-Conventional Transmissible Agents (prions), This cycle must be completed by a period of 60 min disinfection for the instruments concerned. Cf. Recommendations of the instruction DGS/R3/2011/N/449 of 1st December 2011

HANDLING/STORAGE

When handling of sterile boxes, special care must be paid so to avoid damaging the sterile barrier The health care establishment must establish a storage period for packaged instruments, depending on the sterile packaging and the recommendations of the manufacturer of this packaging The health care establishment must store sterile ancillaries in a suitable location which is not close to the sterile barrier to be prevented Instruments and ancillaries must not be exposed to direct sunlight, to ionising radiation or to extreme temperatures. They must be stored in a clean, ventilated room It is essential that a sterile ancillary can be distinguished from a non-sterile ancillary ADSM prohibits the transportation of sterile ancillaries outside the health care facility and disclaims any liability for such transportation.

REFERENCES

29CFR1910.1030: Bloodborne pathogen ISO 17663-1: Sterilisation of health care products - Information to be provided by the manufacturer for the processing of re-sterilizable medical devices AAMI TIR12: Designing, testing, and labelling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers. AAMI TR 30: A compendium of processes, materials, test methods, and acceptance criteria for clearing reusable devices ANSIAAM 5729: Comprehensive guide to steam sterilisation and sterility assurance in health care facilities and in non-clinical settings ISO 17665-1: Sterilisation of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilisation process for medical devices (EN 15883-1: Washer-disinfectors – Part 1: General requirements, terms and definitions and tests.

ISO 15883-2: Washer-disinfectors – Part 2: Requirements and tests for washer-disinfectors employing thermal disinfection for surgical instruments, anaesthetic equipment, bowls, dishes, receivers, utensils, glassware, etc.

ADDITIONAL INFORMATION

For any additional information about the device or to enquire about specific surgical techniques, kindly contact ADSM's client services or the distributor. The instruction for use can be downloaded from the www.adsm-group.com/eli/. If the website is unavailable or you face difficulties to access the document, ADSM can send the instruction for use either by e-mail in PDF format or by post at no extra cost (see contact details above for sending us your request).

PRODUCT COMPLAINTS

Any customer or user of this system who has a complaint or any reason for dissatisfaction (defective quality of a product, a defect, efficacy, durability, reliability, safety, effectiveness and/or its performance not function well, or if poor functioning is suspected, the distributor or ADSM must be notified about this. If one of ADSM's products is ever suspected of having caused a patient's death or serious injury, ADSM or the distributor and the competent authority of the Member State in which the user and/or patient is established should be informed about this immediately by telephone, fax or written correspondence. For any complaint, please indicate the name, reference number and batch number of the components and the circumstances of the complaint, the nature of your complaint, and specify whether a written report is requested from the distributor or ADSM. Any return of instruments to ADSM or a distributor must be accompanied by a certificate guaranteeing that they have been cleaned, disinfected and sterilised beforehand.

INFORMATION TO BE PROVIDED TO THE PATIENT *(only for long-term use devices)*

STERILEBIZ & MISTERI devices for long-term use are designed for fixation and stabilization of bone segments for bone fusion. Specific indications are described in the paragraph: "INDICATIONS, CONTRAINDICATIONS AND POSSIBLE UNDESIRABLE EFFECTS". This device does not restore the functionality that can be expected of healthy bone and the patient should not expect unrealistic functional results. In addition, the anatomy of the human body limits the size of any artificial consolidation device used in surgery. This geometrical fact increases the possibility of the onset of mechanical complications such as dismantling, deformation or rupturing of the device. Any complication may lead to an additional surgical procedure to remove the device or potentially insert another device. Consequently, it is very important that you carefully follow your doctor's postoperative specific advice. It is recommended that you do not limit to those you are given your surgeon's advice. Immobilisation and other devices that partially or fully support your weight must be used on your medical practitioner's advice. By following these instructions, you will increase your chances of achieving the desired result, and you will be reducing the risk of injury and/or infection. The summary of safety and clinical performance of STERILEBIZ & MISTERI Kirschner wires is available on request from the manufacturer (in the absence of operation of the EUADMED database). Basic UDI-DI = 37010910D007K7VSTPH & 37010910D007K7WNONSTELZ

REMOVAL AND DISPOSAL OF MEDICAL DEVICES

The removal and handling of surgical devices is to be done in accordance with the recommendations of the standard ISO 12891-1 "Removal and Analysis of Surgical Implants, Part 1: Removal and Handling". The surgical technique details the surgical steps related to removal of the device. Explanted devices must be disposed of in compliance with the applicable legislation in the country of use.

Special attention will be paid to all the devices in the STERILEBIZ & MISTERI ranges because they have sharp edges. It is strictly forbidden to manipulate them by their extremities. They will have to be placed in specific reprocessing containers avoiding any risk of perforation and thus of injury of the staff.

No specific provision is required for the disposal of unused medical devices. **MAGNETIC RESONANCE (only for long-term use devices)** Magnetic Resonance Imaging (MRI) is a medical imaging technique that provides non-invasive views of the body interior using a strong magnetic field. Consequently, any ferromagnetic metallic foreign body is to be proscribed which can cause burns or injuries by displacement of the body under the effect of attraction of the magnetic field. Stainless steel and Phynox used for devices of the STERILEBIZ & MISTERI ranges are non-ferromagnetic metallic materials which pose no known risk for possible migration due to exposure to the magnetic resonance environment. It is the same for the healing for which the literature found compared to system has been concluded negligible effects of heating also for multiple devices closely implanted. A delay of 6-8 weeks is recommended after implantation before performing this type of examination. However, the system has not been tested for safety and compatibility in the MR environment. The system has not been tested for heating or migration in the MR environment. Thus, it is possible that this product will cause interference with clinical investigation tools such as MRIs. Some artifacts can be observed.







Conditional 8 of ASTM F2503:

Note: This information pertains to a device that has MRI labeling at 1.5-Tesla and 3-Tesla, only. The device was determined to be MR Conditional per ASTM F2503 (Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment). Non-clinical testing demonstrated that the medical device is MR Conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only**
- Maximum spatial gradient magnetic field of 20-gauss/cm (a higher spatial gradient magnetic field value may apply if properly calculated)**
- Maximum MR system specific absorption rate (SAR) of 2-W/kg (i.e., Normal Operating Mode) for 15 minutes of scanning (per pulse sequence)**

WARNINGS

The manufacturer's responsibility is limited to the instructions and applications included in this instruction manual.

SYMBOL	MEANING
	DO NOT RESTERILISE
	USE BY DATE
	BATCH CODE
	STERILISATION BY RADIATION - DOUBLE STERILE BARRIER SYSTEM
	NON STERILE DEVICE
	CAUTION
	MANUFACTURER
	DO NOT USE IF PACKAGING DAMAGED
	KEEP AWAY FROM SUNLIGHT
	KEEP DRY
	CONSULT ELECTRONIC INSTRUCTIONS FOR USE
	REFERENCE
	MATERIAL
	DIAMETER
	CAUTION - FEDERAL, USA LAW RESTRICTS THE SALE AND/OR THE USE OF THIS DEVICE TO OR ON THE ORDER OF A PHYSICIAN
	NUMBER OF UNITS IN THE PACKAGING
	COMPATIBLE WITH MRI UNDER CERTAIN CONDITIONS
	MEDICAL DEVICE
	UNIQUE DEVICE IDENTIFIER
	DATE OF MANUFACTURE + COUNTRY OF ORIGIN
	SWISS AUTHORISED REPRESENTATIVE + IMPORTER
	DISTRIBUTOR
	CE MARK
	UK AUTHORIZED REPRESENTATIVE

UK REP

Ortho Consulting Group
Building 81, Greenham Park Park,
Berkshire, Reading, RG2 9AT,
United Kingdom

ADSM
7b rue Lavoisier, 69680 Chassieux-France
Tel: +33 (0)4 28 71 03 10 / Fax: +33 (0)4 28 71 03 20
office@synchromedical.com

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