to check that they are not damaged. Any medical device showing signs of damage must not be used. Additional information about these medical devices can be obtained from ADSM or the distributor.

PROCESSING (Non-sterile delivered device) or REPROCESSING (reusable device for multiple patients)

It is incumbent to health care facilities to ensure that the re-treatment is done using appropriate equipment and materials and that the personnel responsible for the re-treatment have been correctly trained in order to obtain the desired result. The equipment and the processes must be validated and regularly monitored.

WARNINGS
Surgical instrumentation must be manipulated and/or used by trained, qualified persons having read the present instruction notice

Surgical instrumentation must be used in premises compatible with standard aseptic conditions and with current practices in health care establishments
The ancillaries and instruments are delivered as non-STERILE.
The transfer of an ancillary between two health care establishments is forbidden without prior return to ADSM

Certain ancillaries are constituted by temporary invasive devices such as drills, reamers, drills, taps and probes. These instruments must be used with precaution.
The equipment must only be used for the function for which they are dedicated in operating techniques. Any deterioration of the equipment can lead to risks of malfunction

GENERALITIES
Before being used, the instruments, whether new or not, must be treated carefully according to the instructions
The instruments boxes do not constitute a sterile barrier; a sterile packaging must be used to keep them in a sterile state

Ancillary equipment is delivered non-sterile in transport containers. A transfer sheet (which specifies the non-sterile state of the ancillary) is supplied with each ancillary for checking at reception and then before sterilisation
Instruments or ancillaries generally have a long service life; however, poor handling or unsuitable protection can rapidly reduce it. Instruments which are no longer efficient because of long use, poor handling or inappropriate maintenance, must be returned to ADSM

Any malfunction concerning the instruments must be indicated to ADSM
ADSM instruments and ancillaries are delivered in their original sterile packaging, cleaning, disinfection, inspection and final sterilisation before being sent back to ADSM. The documentation on decontamination must be supplied with instruments returned to ADSM
Any non-sterile medical device received must be inspected upon receipt to verify the integrity and the proper functioning of the device.
Washer-disinfectors used by the health care facility must comply with ISO 15883.

CLEANING

Caution:
The following precautions must be observed by hospital personnel who work with contaminated or potentially contaminated medical devices. Pointed or cutting devices must be handled with great care

Instruments or metal brushes or scouring pads during manual cleaning procedures. These items risk damaging the surface and the finish of the instruments. Use nylon brushes with soft bristles and swabs with dimensions adapted to the devices to treat

Do not allow contaminated devices to dry before re-treatment so as to facilitate later steps
Saline solution and cleaning/disinfecting agents containing aldehyde, mercury, active chlorine, chlorine, bromine, iodine or iodide are corrosive and must not be used
Do not use mineral oil or silicone lubricants as they coat the micro-organisms preventing direct contact of the steam with the surface and are difficult to eliminate.

Cleaning and neutral enzymatic agents are recommended and preferred for cleaning reusable ADSM devices
Only agents with a neutral pH must be used for cleaning aluminium devices. Contact with strong alkaline detergents or with strong oxidising caustic soda, iodine or chlorine must be avoided as aluminium and stainless steel can be chemically attacked and the device damaged. The user must always refer to and comply with the instructions provided by the cleaning product manufacturer

The use of hard water must be avoided. Softened tap water can be used for initial rinsing. Final rinsing must be done with purified water in order to eliminate deposits on the instruments. One or several of the following processes can be employed to purify water: ultrafiltration (UF), reverse osmosis (RO), de-ionisation or equivalent

Comment 1: Drills, broaches, rasps and cutting instruments must be carefully inspected after treatment with alkaline detergents in order to ensure that the cutting edges are not deteriorated.

Comment 2: In the case of first use of the instruments, the establishment shall perform a general cleaning and disinfection procedure in order to ensure the initial condition of the instruments. A minimum of 3 consecutive cleaning operators is recommended.

Instructions

Where the device will be used

Remove excess liquid and organic tissues on the instruments with a lint-free disposable pad. Place the devices in a tray of distilled water and cover with a damp cloth
Remove the plastic film and potential protective foam
Visually check that all the instruments cited on the attached check-list are present
Comment: the instruments must be cleaned within 30 minutes following their use to minimise the possibility of drying before cleaning.

Make sure that the instruments are completely dry before cleaning.
The instruments used must be sent to the procurement service in closed or covered containers in order to avoid any new risk of contamination
When this is applicable, the instruments with several components must be disassembled for cleaning and drying. Make sure that the components are completely dry.
When this is applicable, articulated instruments must be open
All cleaning agents must be prepared observing the dilution and the temperature recommended by the manufacturer. Softened tap water can be used for preparing the cleaning agents.
The use of recommended temperatures is important for optimal performances of the cleaning agents

Comment: new cleaning solutions must be prepared.

Procedure for automatic cleaning/disinfection

Soaking (pre-cleaning)
Completely immerse the instruments in an alkaline solution having fungicide, bactericides and virucide powers and leave them to soak for at least 15 minutes at room temperature. We recommend the use of the product NIODISHERSEPTO PRECLEAN ZP (CJ) = 5 ml / l of the brand DR WEIGERT Use a nylon, flexible bristle brush to softly brush down each unit and all visible stains are eliminated. Pay particular attention to asperities, cavities, contact surfaces, connectors and holes in order to clean areas. Act on the articulated parts in order to get to difficult access areas

Comment: it is not recommended to prolong the soaking step for an indeterminate period at the risk of deteriorating the instruments.
Remove the device from the alkaline solution and rinse abundantly with lukewarm water of the network for at least 5 minutes. Carefully and effectively rinse cavities, holes and any other difficult access areas

Place the instruments in the tray of a suitable washing machine/autoclave and treat using the standard cycle of the instrument washing machine/autoclave as per the instructions from the manufacturer of the washing machine/autoclave
Washing with a suitable alkaline solution suitable for the washing machine/autoclave (observe the concentrations and the temperatures of use of the products). We recommend the use of the product NIODISHERSEPTO CLEAN (CJ)=5mL/ DR WEIGERT's brand.

Thermal disinfection assuring a value of A greater than 3000, i.e. A minimum of 3 min at 93°C (or any other equivalent time/temperature pair)

Rinse with demineralised water
Dry in a clean, dry environment. We recommend the use of MEDIKLAR® SPECIAL drying activator (CJ) = 0.8 ml / l of the DR WEIGERT brand.
Comment: The ancillary equipment is designed to operate at room temperature. It is recommended not to use the equipment directly after steam sterilisation and to wait for the temperature to return to ambient. It is the responsibility of the user, in consultation with the personnel of the operating theatre to check the proper coupling of the ADSM instruments with the equipment available in the operating theatre.

Inspection before packaging
Carefully inspect each device to ensure the elimination of any visible contamination. If there is any contamination, repeat the cleaning/disinfection process
Check the action of movable parts (e.g. hinges, connectors, sliding parts, springs, etc.) over their complete amplitude

Check the absence of deformation of the instruments, in particular instruments which are connected to motors
When the instruments are assembled, check that the devices fit together correctly with any components in contact

Comment that there is no excessive play between the assembled components

PACKAGING

Medical instrument packaging or steam sterilisation sachets available commercially, can be used for wrapping individual instruments. The wrapping must be prepared using a protective double wrapping technique or an equivalent method

For sterilisation, trays and boxes with covers can also be placed in an approved sterilisation container with sealing lid. Follow the sterilisation container manufacturer's instructions for insertion and the replacement of sterilisation filters in sterilisation containers

Instrument transport and storage boxes and trays must be used in the following conditions: All devices must be properly arranged in order to ensure contact of the steam with all instrument surfaces. The instruments must not be stacked or placed in close contact (use the arranging supports or silicone mat provided for this)
The user must check that the contents of the instrument box do not spill out, once the devices are arranged in the box

Only devices manufactured and/or distributed by ADSM must be placed in ADSM trays

STERILIZATION

The cleaning/decontamination step is compulsory before sterilising. It is important that sufficient cleaning is done before sterilisation without which, the sterilisation can be compromised

After cleaning, the instruments must be rapidly and meticulously reassembled and placed in their respective trays/locations in order to avoid any possible damage or contamination
The instrument trays and supports must have been decontaminated and cleaned according to the same arrangements as those for the instruments

The health care establishment is responsible for the internal procedures for reassembly, inspection and packaging of instruments in order to ensure the sterilising penetration of the steam and appropriate drying. The arrangements to make for the protection of the patient, especially injuring edges of instruments must also be recommended by the health care establishment

Sterilisation by steam/moist heat is the recommended method for the ADSM set of surgical instruments. Refer to the table below for the parameters of the sterilisation cycles validated by ADSM ensuring a sterility assurance level (SAL) of 10⁻⁶

Sterilisation methods using ethylene oxide or plasma must not be used
The recommendations of the manufacturer of the steriliser must always be followed. During the sterilisation of several instrument sets in a sterilisation cycle, ensure that the maximum load indicated by the manufacturer is not exceeded
The set of instruments must be correctly prepared and packaged in the trays and/or boxes to allow the steam to penetrate and enter in direct contact with all surfaces

Validated parameters for steam sterilization:					
	Type of cycle	Temperature (CJ) (°F)	Exposure Time	Drying Time	
US	Pre-vacuum	270°F	4 min	20 min	
International ¹⁾	Pre-vacuum	134°C / 273°F	3 min	20 min	
World Health Organization (WHO)	Pre-vacuum	273°F	18 min	20 min	

NB:

* Steam sterilization cycles with longer times than those listed are also acceptable

** This cycle not recommended in the USA

¹⁾ Recommended by the CDC for the reprocessing instruments where there is concern regarding TSE/CJD and prion diseases. The use of an operation with a risk of Non-Conventional Transmissible Agents (prions), this cycle must be completed by a period of total deactivation for the instruments concerned. Cf. Recommendations of the instruction DGS/R3/2011/Nr149 of 1st December 2011

HANDLING/STORAGE

When handling of sterile boxes, special care must be paid so to avoid damaging the sterile barrier.

The health care establishment must establish a storage period for packaged instruments, depending on the state of packaging and the recommendations of the manufacturer of this packaging

The health care establishment must store sterile ancillaries in a suitable location which will allow damage to the sterile barrier to be prevented

Instruments and ancillaries must not be exposed to direct sunlight, to ionising radiation or to extreme temperatures. They must be stored in a clean, ventilated room
The instruments must be stored in their original packaging, in order to ensure the sterility of ADSM any liability for such transportation.

REFERENCES

The instruction DGS/R3/2011/Nr149 of 1st December 2011 concerning the update of the recommendations aimed at reducing the risks of Non-Conventional Transmissible Agents (prions) during invasive operations
20CFR910.1030: Bloodborne pathogen
ISO 17664: Sterilisation of medical devices -Information to be provided by the manufacturer for the processing of re-sterilisable medical devices

AAMI TIR.12: Designing, testing, and labelling reusable medical devices for re-processing in health care facilities: A guide for medical device manufacturers
AAMI TIR.30: A compendium of processes, materials, tests methods, and acceptance criteria for cleaning reusable devices

ANSI/AAMI S179: Comprehensive guide to steam sterilisation and sterility assurance in health care facilities.
ANSI/AAMI ST.35: Safe handling and biological decontamination of reusable medical devices in health care facilities and in non-clinical settings
ISO 17665-1: Sterilisation of health care products - Moist heat - Part 1: Requirements for the development, validation and routine control of a sterilisation process for medical devices
ISO 15883-1: Washer-disinfectors – Part 1: General requirements, terms and definitions and tests
ISO 15883-2: Washer-disinfectors – Part 2: Requirements and tests for washer-disinfectors employing manual disinfection for surgical instruments, anaesthetic equipment, bowls, dishes, receivers, utensils, glassware, etc.

ADDITIONAL INFORMATION

For any additional information about the device or to enquire about specific surgical techniques, kindly contact ADSM's client services or the distributor.

PRODUCT COMPLAINTS

Any customer or user of this system who has a complaint or any reason for dissatisfaction concerning the quality of a product, its identity, durability, reliability, safety, effectiveness and/or its performance must notify the distributor or ADSM. In addition, should any of the implanted components not function well, or if poor functioning is suspected, the distributor or ADSM must be notified about this. If one of ADSM's products is ever suspected of having caused or contributed towards a patient's death or serious injury, ADSM or the distributor and the competent authority of the Member State in which the user and/or patient is established should be informed about this immediately by telephone, fax or written correspondence. For any complaint, please indicate the name, reference number and batch number of the component(s) concerned, your name and address, the nature of your complaint, and specify whether a written report is requested from the distributor or ADSM. Any return to ADSM or ADSM or a distributor must be accompanied by a certificate guaranteeing that they have been cleaned, disinfected and sterilised beforehand.

INFORMATION TO BE PROVIDED TO THE PATIENT (*only for long-term use devices*)

STERILEBITZ & MISTERI devices for long-term use are designed for fixation and stabilization of the body interior using a strong magnetic field. Consequently, any ferromagnetic metallic foreign body is to be proscribed which can cause burns or injuries by displacement of the body under the effect of attraction of the magnetic field. Stainless steels and Phynox used for devices of the STERILEBITZ & MISTERI ranges are non-ferromagnetic metallic materials which pose no risk of possible migration due to exposure to the magnetic resonance environment it is the same for the heating for which the literature found compared to system has been concluded negligible effects of heating also for multiple devices closely implanted. A delay of 6-8 weeks is recommended after implantation before performing the type of examination. However, the system has not been tested for safety and compatibility in the MR environment. The system has not been tested for heating or migration in the MR environment. Thus, it is possible that this product will cause interference with clinical investigation tools such as MRIs. Some artifacts can be observed.

This device does not restore the functionality that can be expected of healthy bone and the zirconium/zirconium oxide does not expect unrealistic functional results. In addition, the anatomy of the human body limits the size of any artificial consolidation device used in surgery. This geometrical limitation increases the possibility of the onset of mechanical complications such as dismantling, deformation or rupturing of the device. Any complication may lead to an additional surgical procedure to remove the device or potentially insert another device.

Consequently, it is very important that you carefully follow your doctor's postoperative instructions. It is recommended that the activities you do be limited to those on your surgeon's advice. In addition, the use of other devices that are not fully supported your weight must be used on your medical practitioner's advice. By following these instructions, you will increase your chances of achieving the desired result, and you will be reducing the risk of injury and/or additional surgical procedures.

Information should be provided to the patient in accordance with the EUDAMED website : <https://ec.europa.eu/eudamed/index>

Basic UDI-UI = 37010910D7007K5T8PHE & 37010910D7007KXWVN0STEUZ

REMOVAL AND DISPOSAL OF MEDICAL DEVICES

The removal and handling of surgical devices is to be done in accordance with the instructions of the standard ISO 12891-1 "Removal and Analysis of Surgical Implants. Part 1: Removal and Handling". The surgical technique details the surgical steps related to removal of the device.

Expanted devices must be disposed of in compliance with the applicable legislation in the country of use.
Special attention will be paid to all the devices in the STERILEBITZ & MISTERI ranges because they have sharp edges. It is strictly forbidden to manipulate them by their extremities. They will have to be placed in specific reprocessing containers avoiding any risk of perforation and loss of injury of staff.
No specific provision is required for the disposal of unused medical devices.

MAGNETIC RESONANCE (*only for long-term use devices*)

Magnetic Resonance Imaging (MRI) is a medical imaging technique that provides non-invasive view of the body interior using a strong magnetic field. Consequently, any ferromagnetic metallic foreign body is to be proscribed which can cause burns or injuries by displacement of the body under the effect of attraction of the magnetic field. Stainless steels and Phynox used for devices of the STERILEBITZ & MISTERI ranges are non-ferromagnetic metallic materials which pose no risk of possible migration due to exposure to the magnetic resonance environment it is the same for the heating for which the literature found compared to system has been concluded negligible effects of heating also for multiple devices closely implanted. A delay of 6-8 weeks is recommended after implantation before performing the type of examination. However, the system has not been tested for safety and compatibility in the MR environment. The system has not been tested for heating or migration in the MR environment. Thus, it is possible that this product will cause interference with clinical investigation tools such as MRIs. Some artifacts can be observed.

Conditional 8 of ASTM F2503:




Note: This information pertains to a device that has MRI labeling at 1.5-Tesla and 3-Tesla, only.

The device was determined to be MR Conditional as per ASTM F2503 (Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment). Non-clinical testing determined that the medical device is MR Conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only*
- Maximum spatial gradient magnetic field of 720-gauss/cm (a higher spatial gradient magnetic field value may apply if properly calculated)*
- Maximum MR system reported whole-body-averaged specific absorption rate (SAR) of 2-W/kg (i.e., the Normal Operating Mode) for 15 minutes of scanning (per pulse sequence)*

WARNINGS

The manufacturer's responsibility is limited to the instructions and applications included in this instruction manual.

SYMBOL	MEANING
	DO NOT RESTERILISE
	VALID UNTIL
	BATCH CODE
	STERILISATION BY RADIATION - DOUBLE STERILE BARRIER SYSTEM
	NON STERILE DEVICE
	CAUTION: SEE INSTRUCTION MANUAL
	MANUFACTURER
	DO NOT USE IF PACKAGING DAMAGED
	STORE OUT OF DIRECT LIGHT
	KEEP DRY
	READ THE INSTRUCTIONS
	REFERENCE
	MATERIAL
	DIAMETER
	CAUTION : FEDERAL (USA) LAW RESTRICTS THE SALE AND/OR THE USE OF THIS DEVICE ON THE ORDER OF A PHYSICIAN
	NUMBER OF UNITS IN THE PACKAGING
	COMPATIBLE WITH MRI UNDER CERTAIN CONDITIONS
	MEDICAL DEVICE
	UNIQUE DEVICE IDENTIFIER
	DATE OF MANUFACTURE + COUNTRY OF ORIGIN
	SWISS AUTHORISED REPRESENTATIVE + IMPORTER
	DISTRIBUTOR
	CE MARK

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