-FR-

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INFORMATION MEDICALE IMPORTANTE  
**TOEGRIP®CLASSIC et TOEGRIP®EVO**

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#### INDICATIONS, CONTRAINDICATIONS AND POSSIBLE SIDE EFFECTS

**Indications :**
TOEGRIP®CLASSIC and TOEGRIP®EVO implants are indicated for the reconstruction of minor bones limited to the interdigital fusion of the toes on the mature skeleton in the following cases :

- Rigidified deformations:
  - Affecting the proximal interphalangeal joint (PIJ),
  - Hammer toe deformities,
  - Claw toe deformities (PIP and DIP joints).
- Hammer toe revision surgeries
- Shortening osteotomies of the proximal phalanx

#### ATTENTION: to be used by or on the order of a licensed physician

The licensed physician must take note of the documents accompanying the device. (IFU and surgical technique). No specific training is required for the understanding and implantation of the device. The medical doctor's qualifications and the reading of the accompanying documents are sufficient

#### Contraindications:

- Abs non-exhaustive list of contraindications is as follows:
- Any sign of generalised or local infection.
  - Pathological obesity.
  - Pregnancy.
  - Any other medical or surgical condition that may compromise the success of surgery with instruments, such as the presence of malignant tumours, or serious congenital anomalies, an increase in sedimentation rates that cannot be attributed to other diseases, an increase in the number of white blood cells or a downward trend in such blood cells.
  - Suspected or known allergy or intolerance to the implant's component materials.
  - Any situation requiring the use of different materials
  - Any case not listed in the indications.
  - Any patient who is not willing to follow the postoperative instructions.
  - Any patient in whom use of the implant may interfere with anatomical organs or some expected physiological function. The contraindications related to these devices are similar to those related to other interphalangeal devices. These interphalangeal instruments have not been designed for, intended or sold for any use other than those indicated.

#### Possible side effects:

- In addition to risks associated with any interphalangeal arthrodesis surgery, the list of possible undesirable side effects is as follows (non-exhaustive list):
- Migration, fracture / rupture, expulsion, dislocation, or unusable movement of one or more of the components
  - Persistent pain and disability
  - Bone loss related to short-circuiting of the constraints.
  - The body's reaction to a foreign object owing to the presence of the implant (e.g. the development of a tumour, auto-immune disease and / or poor healing)
  - Skin incompatibility with components where there is insufficient tissue coverage with pain and abnormal sensations due to the volume of the material.
  - Superficial or deep infection with inflammatory reaction.
  - Sepsis
  - Pseudarthrosis or delay of bone fusion.
  - Recurrence of the deformation, loss of correction
  - Excessive loss of blood during the operation and / or haematoma.
  - Lesions and / or vascular disorders (thrombosis)
  - Peripheral circulatory complication (thrombosis) and / or pulmonary embolism.
  - Inability to resume the activities of normal everyday life.
  - Compression of nearby tissues or organs.

**NOTE:** Some additional surgical procedure may be required to remedy some of these possible undesirable events.

#### INSERTION OF THE DEVICE

The TOEGRIP®CLASSIC and TOEGRIP®EVO implants must be used with instruments specially designed for the device according to the surgical technique concerned. It must not be used with any other instruments, unless this has specifically been recommended in some other ADSM documentation, as using the system with other instruments poses the risk of incompatibility and cannot be guaranteed. ADSM shall not be held responsible in the case of using instruments that are not related to the implanted system.
**NEVER REUSE AN IMPLANT.**
Reuse of an implant presents the risk of contamination and loss of mechanical integrity.

#### PACKAGING

The packaging of the device must be intact on receipt. Any device that has been damaged or any device whose packaging is damaged must not be used and must rather be returned to ADSM.

#### STORAGE AND HANDLING CONDITIONS

The implants must be stored in their original packaging, out of direct sunlight, humidity and extreme temperatures.

#### MAGNETIC RESONANCE

Magnetic Resonance Imaging (MRI) is a medical imaging technique that provides non-invasive views of the body interior using a strong magnetic field. Consequently, any ferromagnetic metallic foreign body is to be proscribed which can cause burns or injuries by displacement of the body under the effect of attraction of the magnetic field. The PEEK used to manufacture the TOEGRIP®CLASSIC and TOEGRIP®EVO systems is therefore a non-ferromagnetic and non-conductive polymer (plastic). Consequently, it poses no known risk of possible heating or migration due to exposure to the magnetic resonance environment.

#### VERIFICATION

The devices must always be checked before use. Those that present signs of damage or scratches on the surface must not be used. The TOEGRIP®CLASSIC and TOEGRIP®EVO implants are delivered sterile and have been sterilised through gamma radiation.
You need to carefully check that the round disc (sterilisation mark) on the packaging is red. Do not use an implant whose disc isn't red.
You also need to check that the expiry date on the packaging has not passed. ADSM disengages itself from any liability in the case of the use of implants past their expiry date. Any implant with damaged or torn packaging must not be restérilised but must rather be returned to ADSM.
**TAKE NOTE:** Do not restérilise an implant that has already been sterilised.

#### WARNINGS AND PRECAUTIONS

Every surgery carries the risk of failure, as many external factors can compromise the outcome, especially if the patient's state of health is already compromised. The implants must only be used in cases where bone fusion is deemed to be possible. No implant has been designed to withstand body weight in the long term without being accompanied by solid bone fusion. If a delay or lack of consolidation occurs, the implant may migrate and/or rupture. The functional product lifetime is of 6 months. In fact, the first bone fusion occurs around 90 days of keeping proxy. Therefore, from 6 months the product is still implanted but it has no functional performance.

Given the current state of knowledge and device anteriority, the risk caused by aging cannot be totally ruled out. The following of preoperative and operative procedures, knowledge of surgical techniques, the selection of implants and their correct placement as well as an appropriate selection and the patient's cooperation are factors leading to the operation's success. It has been shown that patients who smoke present a high risk of failure to obtain bone fusion. They need to be informed of this and any possible consequences. Obese patients, those who are malnourished and/or consume alcohol as well as patients with poor muscle or bone quality and/or those suffering from paralysis are also poor candidates for interphalangeal osteosynthesis surgery.

-EN-

## INSTRUCTION FOR USE

IMPORTANT MEDICAL INFORMATION  
**TOEGRIP®CLASSIC and TOEGRIP®EVO**

Reference : SUP\_7.015 version 12  
Date of last update : 01/2024  
Date of first CE marking :  
- TOEGRIP®CLASSIC : 2010  
- TOEGRIP®EVO : 2016

Product manufactured by :

 <b>ADSM</b>	 <b>UK REP</b>	 <b>MDS5-UK RP Ltd.</b>
<b>7b rue Lavoisier, 69680 Chassieu – France</b>	<b>0197 2571</b>	<b>0197 2571</b>
<b>Tel<span> </span>: +33 (0)4 28 71 03 10 // Fax<span> </span>: +33 (0)4 28 71 03 20</b>		<b>office@synchromedical.com</b>

**CAUTION: USA Federal Law restricts this device to sale by or on the order of physician.**

#### FDA Specific - CAUTION

#### Specific IFU is dedicated to USA.

Please refer to our website, section "IFU": [www.synchromedical.com](http://www.synchromedical.com).
In case you prefer receiving IFU in a paper-based version, please send an email to [office@synchromedical.com](mailto:office@synchromedical.com) or call +33 428 710 320. IFU will be sent to the address you will provide without any additional fee.

#### PURPOSE

The TOEGRIP®CLASSIC and TOEGRIP®EVO interphalangeal implants are designed to relieve pain and disability of the forefoot by restoring and/or maintaining the alignment of two adjacent digital bone segments to optimize the achievement of a correct bone fusion for the concerned segments. These devices are only intended to be used for the forefoot of a mature skeleton.

#### DESCRIPTION:

The TOEGRIP®CLASSIC implant consists of a single part implant with three flexible intramedullary rods inserted into each part of the phalanx. The rods are attached by a T-shaped structure. The attachment concept is based on a press-fit contact due to its tapered shape and the macrostructures solidly anchored once impacted, thus avoiding any displacement of the implant. The TOEGRIP®CLASSIC implant is intended for single use only and is available in a range of 5 sizes with 3 possible degrees: 0°, 10° or 20°.

The TOEGRIP®EVO implant consists of a single part implant with three flexible intramedullary rods inserted into each part of the phalanx. The rods are connected by a cylindrical structure with a sagittal height offset between the proximal and distal portion. The attachment concept is based on a press-fit contact due to its tapered shape and the macrostructures solidly anchored once impacted, thus avoiding any displacement of the implant. The TOEGRIP®EVO implant is intended for single use only and is available in a range of 5 sizes with 3 possible degrees: 0°, 10° or 20°.

TOEGRIP®CLASSIC and TOEGRIP®EVO are made from PEEK according to Standard ASTM F2026.

#### INDICATIONS, CONTRAINDICATIONS AND POSSIBLE SIDE EFFECTS

**Indications :**
TOEGRIP®CLASSIC and TOEGRIP®EVO implants are indicated for the reconstruction of minor bones limited to the interdigital fusion of the toes on the mature skeleton in the following cases :

- Rigidified deformations:
  - Affecting the proximal interphalangeal joint (PIJ),
  - Hammer toe deformities,
  - Claw toe deformities (PIP and DIP joints).
- Hammer toe revision surgeries
- Shortening osteotomies of the proximal phalanx

#### ATTENTION: to be used by or on the order of a licensed physician

The licensed physician must take note of the documents accompanying the device. (IFU and surgical technique). No specific training is required for the understanding and implantation of the device. The medical doctor's qualifications and the reading of the accompanying documents are sufficient

#### Contraindications:

- Abs non-exhaustive list of contraindications is as follows:
- Any sign of generalised or local infection.
  - Pathological obesity.
  - Pregnancy.
  - Any other medical or surgical condition that may compromise the success of surgery with instruments, such as the presence of malignant tumours, or serious congenital anomalies, an increase in sedimentation rates that cannot be attributed to other diseases, an increase in the number of white blood cells or a downward trend in such blood cells.
  - Suspected or known allergy or intolerance to the implant's component materials.
  - Any situation requiring the use of different materials
  - Any case not listed in the indications.
  - Any patient who is not willing to follow the postoperative instructions.
  - Any patient in whom use of the implant may interfere with anatomical organs or some expected physiological function. The contraindications related to these devices are similar to those related to other interphalangeal devices. These interphalangeal instruments have not been designed for, intended or sold for any use other than those indicated.

#### Possible side effects:

- In addition to risks associated with any interphalangeal arthrodesis surgery, the list of possible undesirable side effects is as follows (non-exhaustive list):
- Migration, fracture / rupture, expulsion, dislocation, or unusable movement of one or more of the components
  - Persistent pain and disability
  - Bone loss related to short-circuiting of the constraints.
  - The body's reaction to a foreign object owing to the presence of the implant (e.g. the development of a tumour, auto-immune disease and / or poor healing)
  - Skin incompatibility with components where there is insufficient tissue coverage with pain and abnormal sensations due to the volume of the material.
  - Superficial or deep infection with inflammatory reaction.
  - Sepsis
  - Pseudarthrosis or delay of bone fusion.
  - Recurrence of the deformation, loss of correction
  - Excessive loss of blood during the operation and / or haematoma.
  - Lesions and / or vascular disorders (thrombosis)
  - Peripheral circulatory complication (thrombosis) and / or pulmonary embolism.
  - Inability to resume the activities of normal everyday life.
  - Compression of nearby tissues or organs.

**NOTE:** Some additional surgical procedure may be required to remedy some of these possible undesirable events.

#### INSERTION OF THE DEVICE

The TOEGRIP®CLASSIC and TOEGRIP®EVO implants must be used with instruments specially designed for the device according to the surgical technique concerned. It must not be used with any other instruments, unless this has specifically been recommended in some other ADSM documentation, as using the system with other instruments poses the risk of incompatibility and cannot be guaranteed. ADSM shall not be held responsible in the case of using instruments that are not related to the implanted system.
**NEVER REUSE AN IMPLANT.**
Reuse of an implant presents the risk of contamination and loss of mechanical integrity.

#### PACKAGING

The packaging of the device must be intact on receipt. Any device that has been damaged or any device whose packaging is damaged must not be used and must rather be returned to ADSM.

#### STORAGE AND HANDLING CONDITIONS

The implants must be stored in their original packaging, out of direct sunlight, humidity and extreme temperatures.

#### MAGNETIC RESONANCE

Magnetic Resonance Imaging (MRI) is a medical imaging technique that provides non-invasive views of the body interior using a strong magnetic field. Consequently, any ferromagnetic metallic foreign body is to be proscribed which can cause burns or injuries by displacement of the body under the effect of attraction of the magnetic field. The PEEK used to manufacture the TOEGRIP®CLASSIC and TOEGRIP®EVO systems is therefore a non-ferromagnetic and non-conductive polymer (plastic). Consequently, it poses no known risk of possible heating or migration due to exposure to the magnetic resonance environment.

#### VERIFICATION

The devices must always be checked before use. Those that present signs of damage or scratches on the surface must not be used. The TOEGRIP®CLASSIC and TOEGRIP®EVO implants are delivered sterile and have been sterilised through gamma radiation.
You need to carefully check that the round disc (sterilisation mark) on the packaging is red. Do not use an implant whose disc isn't red.
You also need to check that the expiry date on the packaging has not passed. ADSM disengages itself from any liability in the case of the use of implants past their expiry date. Any implant with damaged or torn packaging must not be restérilised but must rather be returned to ADSM.
**TAKE NOTE:** Do not restérilise an implant that has already been sterilised.

#### WARNINGS AND PRECAUTIONS

Every surgery carries the risk of failure, as many external factors can compromise the outcome, especially if the patient's state of health is already compromised. The implants must only be used in cases where bone fusion is deemed to be possible. No implant has been designed to withstand body weight in the long term without being accompanied by solid bone fusion. If a delay or lack of consolidation occurs, the implant may migrate and/or rupture.

The functional product lifetime is of 6 months. In fact, the first bone fusion occurs around 90 days of keeping proxy. Therefore, from 6 months the product is still implanted but it has no functional performance.

Given the current state of knowledge and device anteriority, the risk caused by aging cannot be totally ruled out. The following of preoperative and operative procedures, knowledge of surgical techniques, the selection of implants and their correct placement as well as an appropriate selection and the patient's cooperation are factors leading to the operation's success. It has been shown that patients who smoke present a high risk of failure to obtain bone fusion. They need to be informed of this and any possible consequences. Obese patients, those who are malnourished and/or consume alcohol as well as patients with poor muscle or bone quality and/or those suffering from paralysis are also poor candidates for interphalangeal osteosynthesis surgery.

#### OTHER PREOPERATIVE, PEROPERATIVE AND POSTOPERATIVE WARNINGS

#### Implant selection

Selecting the suitable size of the implant for each patient is a crucial factor in the success of the operation, and it is important to refer to the surgical technique in this regard. Once implanted, the implant is subjected to repeated strain, and its strength is limited by adaptation of its geometrical structure to the size and shape

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