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INFORMATION MEDICALE IMPORTANTE TIGERTRACK® system

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- IMPORTANT MEDICAL INFORMATION
TIGERTRACK® system**
- Reference: SUP_7.017 version 3
Date of last update: 03/2019
- Date of first EC marking: 2010
- Product manufactured by:
-
- 6 rue des Merisiers, 68920 Wettolsheim les Erlen – France
Tel: +33 (0)3 89 23 75 64 // Fax: +33 (0)9 70 19 01 32
adv@synchromedical.com
- CAUTION: USA Federal law restricts this device to sale by or on the order of physician.**

INSTRUCTION FOR USE

- IMPORTANT MEDICAL INFORMATION
TIGERTRACK® system**
- Reference: SUP_7.017 version 3
Date of last update: 03/2019
- Date of first EC marking: 2010
- Product manufactured by:
-
- 6 rue des Merisiers, 68920 Wettolsheim les Erlen – France
Tel: +33 (0)3 89 23 75 64 // Fax: +33 (0)9 70 19 01 32
adv@synchromedical.com
- CAUTION: USA Federal law restricts this device to sale by or on the order of physician.**

- INTENDED USE**
The TIGERTRACK® screws system is intended for small bone extremity and largebone bone reconstruction fixation/arthrodesis procedures for general use in skeletally mature individuals. The TIGERTRACK® screws system is intended to be permanently implanted without any other additional device and are delivered sterile.
- PURPOSE**
TIGERTRACK® screws are designed to attach and stabilise bone segments in elective osteotomies on metatarsal bones and phalanges of the foot, in order to optimise the correct bone fusion of the segments concerned. These devices are only intended to be used in the forefoot.
- DESCRIPTION:**
The TIGERTRACK® screw system is intended for single use and is composed of TIGERTRACK® self-compressive headless screws (TIGER22, TIGER29 and TIGER45), TIGERTRACK® snap-off screws (TIGER18 and TIGER25) and Weil's snap-off screws TIGERW.
- TIGERTRACK® self-compressive headless screws (TIGER22, TIGER29 and TIGER45) comprise a single part implant. These self-compressive headless screws are available in several different diameters and lengths.
- Snap-off TIGERTRACK® screws (TIGER18 and TIGER25) comprise a single part implant fitted with a threaded portion and a snap-off connector portion with a wider diameter. These are self-compressive, headless screws.
- Weil's snap-off TIGERTRACK® screws (TIGERW) comprise a single part implant fitted with a threaded portion and a snap-off connector portion with a wider diameter. These are self-compressive screws with a head.

- TIGERTRACK® screws are manufactured in **Ti-6Al-4V ELI titanium according to Standard ISO 5832-3 / ASTM F 136**. With the exception of the TIGERW which is manufactured in a chrome-cobalt alloy according to ISO 5832-7.
- INDICATIONS, CONTRAINDICATIONS AND POSSIBLE SIDE EFFECTS**
Indications:
The TIGERTRACK® screw system is indicated as an attachment system for small bone fractures or minor bone reconstruction in the following cases:
- monocortical or bicortical osteotomy of the foot or hand
- metatarsal or distal or proximal metatarsal osteotomy
- fusion of the first metatarsophalangeal and interphalangeal joint - attachment in the case of treatment of Hallux Valgus (e.g. scarf osteotomy, chevron osteotomy, etc.)
- Akin osteotomy
- Weil osteotomy
- CAUTION: to be used by or on the order of a licensed physician.**
- TAKE NOTE:** To be used on or following the recommendations of a medical practitioner. This device may only be used by a doctor who has in-depth knowledge of the device, its intended use, and related instruments and surgical techniques.

- Contraindications:**
A non-exhaustive list of contraindications is as follows:
1. Any sign of generalised or local infection.
2. Pathological obesity.
3. Pregnancy.
4. Any other medical or surgical condition that may compromise the success of surgery with instruments, such as the presence of malignant tumours, or serious congenital anomalies, an increase in sedimentation rates that cannot be attributed to other diseases, an increase in the number of white blood cells or a downward trend in such blood cells.
5. Suspected or known allergy or intolerance to the implant's component materials.
6. Any situation requiring the use of different materials
7. Any case not listed in the indications.
8. Any patient who is not willing to follow the postoperative instructions.
9. Any patient in whom use of the implant may interfere with anatomical organs or some expected physiological function.
10. Any dental treatment that could affect the passivation layer of metal implants (such as fluoride treatment)

- The contraindications pertaining to these devices are similar to those pertaining to other osteosynthesis instruments. These instruments have not been designed for, intended or sold for any use other than those indicated.
- Possible side effects:**
In addition to risks associated with any surgery to the extremities, the list of possible undesirable side effects is as follows (non-exhaustive list):
1. Migration, fracture/rupture, expulsion, dislocation, or undesirable movement of one or more of the components.
2. Persistent pain and disability.
3. Bone loss related to short-circuiting of the constraints.
4. The body's reaction to a foreign object owing to the presence of the implant (e.g. the development of a tumour, auto-immune disease and/or poor healing).
5. Skin incompatibility with components where there is insufficient tissue coverage with pain and abnormal sensations due to the volume of the material.
6. Superficial or deep infection with inflammatory reaction.
7. Sepsis
8. Pseudarthrosis or delay of bone fusion.
9. Recurrence of the déformation, loss of correction.
10. Excessive loss of blood during the operation and/or haematoma.
11. Lesions and/or vascular disorders (thrombosis)
12. Peripheral circulatory complication (thrombosis) and/or pulmonary embolism.
13. Inability to resume the activities of normal everyday life.
14. Compression of nearby tissues or organs.

- WARNING:** Some additional surgical procedure may be required to remedy some of these possible undesirable events.

- INSERTION OF THE DEVICE**
The TIGERTRACK® screw system must be used with instruments specially designed for the device according to the surgical technique concerned. It must not be used with any other instruments, unless this has specifically been recommended in some other ADSM documentation, as using the system with other instruments poses the risk of incompatibility, and cannot be guaranteed. ADSM shall not be held responsible in the case of using instruments that are not associated with the implanted system.
- NEVER REUSE AN IMPLANT.**
Reuse of an implant presents the risk of contamination and loss of mechanical integrity.

- PACKAGING**
The packaging of the device must be intact on receipt. Any device that has been damaged or any device whose packaging is damaged must not be used, and must rather be returned to ADSM.
- STORAGE AND MAINTENANCE CONDITIONS**
The implants must be stored in their original packaging, out of direct sunlight, humidity and extreme temperatures.
- MAGNETIC RESONANCE**
The TIGERTRACK® screw system has not been evaluated for safety and compatibility in the MR environment. The TIGERTRACK® screw system has not been tested for heating or migration in the MR environment.

- CHECKING**
The devices must always be checked before use. Those that present signs of damage or scratches on the surface must not be used. TIGERTRACK® implants are delivered sterile, and have been sterilised through gamma radiation. You need to carefully check that the round disc (sterilisation mark) on the packaging is red. Do not use an implant whose disc isn't red. You also need to check that the expiry date on the packaging has not passed. ADSM disengages itself from any liability in the case of the use of implants past their expiry date. Any implant with damaged or torn packaging must not be resterilised and must rather be returned to ADSM. TAKE NOTE: Do not resterilise an implant that has already been sterilised.

- WARNINGS AND PRECAUTIONS**
Every surgery carries the risk of failure, as many external factors can compromise the outcome, especially if the patient's state of health is already compromised. The implants must only be used in cases where bone fusion is deemed to be possible. No implant has been designed to withstand body weight in the long term without being accompanied by solid bone fusion. If a delay or lack of consolidation occurs, the implant may migrate and/or rupture. The following of preoperative and operative procedures, knowledge of surgical techniques, the selection of implants and their correct placement as well as an appropriate selection and the patient's cooperation are factors leading to the operation's success. It has been shown that patients who are obese present a high risk of failure to obtain bone fusion. They need to be informed of this and any possible consequences. Obese patients, those who are malnourished and/or consume alcohol as well as patients with poor muscle or bone quality and/or those suffering from paralysis are also poor candidates for osteosynthesis surgery.

- OTHER PREOPERATIVE, PEROPERATIVE AND POSTOPERATIVE WARNINGS**
Implant selection
Selecting the suitable size of the implant for each patient is a crucial factor in the success of the operation, and it is important to refer to the surgical technique in this regard. Once implanted, the implant is subjected to repeated strain, and its strength is limited by adaptation of its geometrical structure to the size and shape of human bone. In order to minimise this strain and to avoid compromising the desired bone fusion, it is important to pay careful attention to the patient's selection criteria, the correct placement of the implant, and postoperative care. In fact, conversely, these constraints can lead to excessive strain placed on the material, resulting in deformation, rupture or loosening of the device likely to cause damage or the need to remove the implant prematurely. Any use in body areas not recommended by ADSM cannot be guaranteed.

- Preoperative precautions**
1. Only patients who meet the criteria described in the indications should be selected.
2. Patients who meet the criteria described in the contraindications should not be selected.
3. The implants must be handled and stored with the greatest care. They should not be scratched or damaged.
4. Before carrying out the procedure, the surgeon must check availability of the range of implants and associated instruments. They must personally handle all the devices before use in order to familiarise themselves with insertion of the implant.
5. Additional sterile components must be available in the event of unexpected need.

- Preoperative precautions**
1. Instruction concerning the surgical technique should be followed closely. Any breakage or improper use of the implant may injure the patient or surgical team.
3. Specific instrumentation must be used to insert the device.
4. Take note: Do not use a device of the incorrect size (diameter and/or length), as this may damage the surrounding tissue.
6. Where possible and necessary, an imaging system may be used to check on placement of the device.

- Preoperative precautions**
The instructions and postoperative warnings provided to the patient by the doctor, and the patient's ability to follow them are extremely important, so as not to compromise the desired bone fusion.
1. Detailed instructions on the device's limits need to be provided to the patient. The patient must be informed that any deformity, loosening, and/or rupture of the device are complications that may occur during the postoperative rehabilitation phase following early use or excessive muscle activity. The risk of the onset of these complications could be heightened if the patient is active, frantically ill or incapable of using some immobilisation device or external support. The patient must be advised in order to prevent falls, jumps or hard impacts on the area operated on or about any other action that could compromise the desired bone fusion.
2. To increase the likelihood of success of the surgery, the patient must not expose themselves to excessive strain which could cause weakening or rupture of the device. The patient must be informed of this risk, and be advised so they are aware of how to limit their physical activity, in the case of rupture, the device will need to be adjusted and/or removed immediately, before any serious injury occurs.
3. The patient should be informed about avoiding smoking or consuming alcohol during the bone consolidation process.
4. Lack of bone fusion would lead to excessive, repeated strain placed on the implant. Through the mechanism of fatigue, this strain could end up causing some deformation, loosening or rupturing of the device. It is important to immobilise the fusion area and to check consolidation through x-ray. Should lack of bone fusion persist, or in the case of rupture, the device will need to be adjusted and/or removed immediately, before any serious injury occurs.
5. Any implant collected following ablation must be treated in a manner to prevent its reuse in some other surgical procedure, in accordance with the regulations described in the removal and disposal section for medical devices contained in this manual.

- ADDITIONAL INFORMATION**
For any additional information about the device or to enquire about specific surgical techniques, kindly contact ADSM's client services or the distributor.
- PRODUCT COMPLAINTS**
Any customer or user of this system who has a complaint or any issue for dissatisfaction concerning the quality of a product, its identity, durability, reliability, safety, effectiveness and/or its performance must notify the distributor or ADSM. In addition, should any of the implanted components not function well, or if poor functioning is suspected, the distributor or ADSM must be notified about this. If one of ADSM's products is ever suspected of having caused or contributed towards a patient's death or serious injury, ADSM or the distributor should be informed of this immediately by telephone, fax or written correspondence. For any complaint, please indicate the name, reference number and batch number of the component(s) concerned, your name and address, the nature of your complaint, and specify whether a written report is requested from the distributor or ADSM.

- INFORMATION TO BE PROVIDED TO THE PATIENT**
The TIGERTRACK® screw system is designed for attachment and stabilisation of bone segments during elective osteotomies of metatarsal bones and phalanges of the foot. The device has only been designed for use in the forefoot. Specific instructions for the implant are described in the paragraph: "INDICATIONS, CONTRAINDICATIONS AND POSSIBLE UNDESIRABLE EFFECTS". This device does not restore the functionality that can be expected of healthy bone and the patient should not expect unrealistic functional results. In addition, the anatomy of the human body limits the size of any artificial consolidation device used in surgery. This geometrical limitation increases the possibility of the onset of mechanical complications such as dismantling, deformation or rupturing of the device. Any complication may lead to an additional surgical procedure to remove the device or potentially insert another device. Consequently, it is very important that you carefully follow your doctor's postoperative instructions. It is recommended that the activities you do be limited to those on your surgeon's advice. Immobilisation and other devices that partially or fully support your weight must be used on your medical practitioner's advice. By following these instructions, you will increase your chances of achieving the desired result, and you will be reducing the risk of injury and/or additional surgical procedures.

- REMOVAL AND DISPOSAL OF MEDICAL DEVICES**
The removal and handling of surgical instruments is to be done in accordance with the recommendations standard ISO 12891-1 "Removal and Analysis of Surgical Implants. Part 1: Removal and Handling". Waste must be disposed of in compliance with the applicable legislation in the country of use. In France, waste is to be handled in accordance with Decree no. 97-1048 of 6 November 1997 relating to the disposal of infectious and related health care waste and anatomical parts. No specific provision is required for the disposal of unused medical devices.

- WARNINGS**
The manufacturer's responsibility is limited to the instructions and applications indicated in this instruction manual.
- | SYMBOL | MEANING |
|--------|--|
| | SINGLE USE |
| | DO NOT RESTERILISE |
| | VALID UNTIL |
| | BATCH CODE |
| | STERILISATION BY RADIATION |
| | CAUTION: SEE INSTRUCTION MANUAL |
| | MANUFACTURER |
| | DO NOT USE IF PACKAGING DAMAGED |
| | STORE OUT OF DIRECT LIGHT |
| | KEEP DRY |
| | READ THE INSTRUCTIONS |
| | REFERENCE |
| | MAT. MATERIAL |
| | Lg. LENGTH |
| | ∅ DIAMETER |
| | CAUTION : FEDERAL (USA) LAW RESTRICTS THE SALE AND/OR THE USE OF THIS ORDER OF A PHYSICIAN |

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