

EXray® 20 PEEK

<u>Implant</u>

PW2060....

Snap-off threaded pin

DESCRIPTION

EXRAY® 20 snap-off threaded pin has been developed as a resectable peek pin. The seat of the implant is prepared with a 2 mm drill. The pin consists of a trocar tip, a threaded portion and a portion with a narrower diameter that can be connected to a power tool.

The EXRAY® 20 has a non-compressive Ø2.5mm continuous thread. The 3 markings are positioned at 15 mm, 20 mm and 25 mm. The full length of the thread is 30 mm. Once implanted, the excess length is removed with a saw blade.

 $\mathsf{EXRAY}^{\$}$ 20 snap-off threaded pins are manufactured in PEEK as per standard ASTM F2026.

- Simple preparation
- Variable screw length
- Used to support stabilization of osteotomies
- Simplified revision
- PEEK threaded pins are radiolucent
- The elasticity modulus is close to the one of bone



Instruments



DRDERING INFORMATION

...L 60 mm

PROTECTIVE BARREL









C € 0197







